AIML - CAPSTONE PROJECT

PNEUMONIA DETECTION CHALLENGE

Mentor: Srihari Nalabolu

Team

NAME	EMAIL	PHONE NUMBER			
Guru Siddarth J	siddarth.mine@gmail.com	+91 9743338960			
Keyur Singh	krajput@ufl.edu	+1 7866816931			
Michael Prasan S	michaelprasansm@gmail.com	+91 9788807095			
M V Pavan kumar	pavanmv.147@gmail.com	+91 9743817778			
Swastika Samanta	sstika1920@gmail.com	+91 9049630396			
Vignesh Jairam	vigirasi@gmail.com	+91 6383442844			

Repository	4
Abstract	4
Problem Statement	4
Objectives	6
Use cases	6
Dataset	7
Implementation and EDA: Importing Packages necessary for the project	8
EDA on stage_2_train_labels File	9
EDA on stage_2_detailed_class_info File	12
Merging of both the dataframes	15
Function to show images from train folder with bounding boxes	24
Showing images without pneumonia	25
Showing images with pneumonia	26
EDA on Metadata of DICOM Images	27
EDA on Modality	29
EDA on Body Part Examined	29
EDA on Understanding Different Positions	29
EDA on Conversion Type	30
EDA on Gender	30
Function to plot graphs on different attributes of the dataframe	30
EDA on PatientSex and Count of patientId - Hue by Class	31
EDA on PatientSex and Count of patientId. Split by Target	31
EDA on PatientAge	31
Distribution on Count of Patients based on PatientAgeBins	32
Image Classification	33
Image Pre-processing	37

Functions for Re-usability:	37
VGG-16 Model	38
Vgg model building	38
ResNet50 Model	41
Resnet50 Model building	42
DenseNet121 Model	44
DenseNet121 Model Building:	46
Object Detection:	49
Object Detection Model Building:	50

Repository

The below link has all the artifacts of the Projects. Following are the documents that are in the repository.

Repo Link: Capstonse-Project-CV1

Artifacts:

- 1. Code of the project in both HTML and IPYNB Formats.
- 2. Submission file in CSV file on the Object Detection Outputs.
- 3. Presentation of the Project Outline.
- 4. Detailed Project Report.

Abstract


- Pneumonia is an infection of one or both of the lungs caused by bacteria, viruses, or fungi.
- There are more than 30 different causes of pneumonia, and they're grouped by the cause. The main types of pneumonia are bacterial, viral, and mycoplasma pneumonia.
- A cough that produces green, yellow, or bloody mucus is the most common symptom of pneumonia. Other symptoms include fever, shaking chills, shortness of breath, low energy, and extreme tiredness.
- Pneumonia can often be diagnosed with a thorough history and physical exam. Tests used to look at the lungs, blood tests, and tests done on the sputum you cough up may also be used.
- Treatment depends on the type of pneumonia you have. Antibiotics are used for bacterial
 pneumonia. It may also speed recovery from mycoplasma pneumonia and some special
 cases. Most viral pneumonias don't have a specific treatment and just get better on their
 own. Other treatments may include a healthy diet, more fluids, rest, oxygen therapy, and
 medicine for pain, cough, and fever control.
- Most people with pneumonia respond well to treatment, but pneumonia can cause serious lung and infection problems. It can even be deadly.
- The aim of this project is to develop a software system to detect the disease Pneumonia using Chest x-ray images.
- This is achieved by using multiple convolutional neural network layers where the chest x-ray images are tested, trained and validated.

Problem Statement

Pneumonia is an infection in one or both lungs. Bacteria, viruses, and fungi cause it. The infection causes inflammation in the air sacs in your lungs, which are called alveoli. Pneumonia accounts for over 15% of all deaths of children under 5 years old internationally. In 2017, 920,000 children under the age of 5 died from the disease. It requires review of a chest radiograph (CXR) by highly trained specialists and confirmation through clinical history, vital

signs and laboratory exams. Pneumonia usually manifests as an area or areas of increased opacity on CXR. However, the diagnosis of pneumonia on CXR is complicated because of a number of other conditions in the lungs such as fluid overload (pulmonary edema), bleeding, volume loss (atelectasis or collapse), lung cancer, or postradiation or surgical changes. Outside of the lungs, fluid in the pleural space (pleural effusion) also appears as increased opacity on CXR. When available, comparison of CXRs of the patient taken at different time points and correlation with clinical symptoms and history are helpful in making the diagnosis. CXRs are the most commonly performed diagnostic imaging study. A number of factors such as positioning of the patient and depth of inspiration can alter the appearance of the CXR, complicating interpretation further. In addition, clinicians are faced with reading high volumes of images every shift.

We aim to achieve this by proposing a model using the Deep Learning method. The dataset for the project is collected from kaggle which contains images of infected and normal lungs. Dataset is divided into 3 categories such as train, test and validate.



Objectives

- To classify patients' chest X-ray images as a pneumonia case or a normal case.
- To build a pneumonia detection system, to locate the position of inflammation in an image
- Assist physicians to make better clinical decisions or even replace human judgment in certain functional areas of healthcare (eg, radiology).
- Guided by relevant clinical questions, powerful AI techniques can unlock clinically relevant information hidden in the massive amount of data, which in turn can assist clinical decision making.

Use cases

- Pneumonia is the single largest infectious cause of death in children worldwide.
 Pneumonia killed 740,180 children under the age of 5 in 2019, accounting for 14% of all deaths of children under five years old but 22% of all deaths in children aged 1 to 5.
- The infection is caused by viruses, bacteria, or fungi and can be prevented by immunization, adequate nutrition, and by addressing environmental factors.
- Pneumonia caused by bacteria can be treated with antibiotics, but only one third
 of children with pneumonia receive the antibiotics they need.
- Detecting pneumonia is a difficult challenge.
- Currently, chest X-rays are one of the best methods for the detection of pneumonia. X-rays are the most common and widely available diagnostic imaging technique, playing a crucial role in clinical care and epidemiological studies.
- It requires review of a chest radiograph (CXR) by highly trained specialists and confirmation through clinical history, vital signs and laboratory exams.
- Pneumonia usually manifests as an area or areas of increased opacity on CXR.
 However, the diagnosis of pneumonia on CXR is complicated because of a
 number of other conditions in the lungs such as fluid overload (pulmonary
 edema), bleeding, volume loss (atelectasis or collapse), lung cancer, or
 post-radiation or surgical changes.
- Outside of the lungs, fluid in the pleural space (pleural effusion) also appears as increased opacity on CXR.

- When available, comparison of CXRs of the patient taken at different time points and correlation with clinical symptoms and history are helpful in making the diagnosis.
- CXRs are the most commonly performed diagnostic imaging study. A number of factors such as positioning of the patient and depth of inspiration can alter the appearance of the CXR, complicating interpretation further.
- In addition to this, experts are faced with the prospect of reading large volumes of images which may prove to be impossible and human error may creep in.
- In 2017,researchers at Stanford came up with a deep learning algorithm capable of diagnosing pneumonia from chest X-ray images. In just over a month of development, the algorithm outperformed expert radiologists in detection.
- Artificial Intelligence can help radiologists and physicians in decision making and can replace human intervention in some cases.
- Algorithms can recognise complex patterns and provide insights and interpretation to big data.
- While machines are better at handling large data and computation, humans can
 be the perfect complement to that by adding the conscious element and use the
 inputs given by the machine for better clinical decisions and treatment.

Dataset

- We will be using the RSNA Pneumonia Detection Challenge for this project.
- This is a two-stage challenge. You will need the images for the current stage provided as stage_2_train_images.zip and stage_2_test_images.zip. You will also need the training data stage_2_train_labels.csv and the sample submission stage_2_sample_submission.csv, which provides the IDs for the test set, as well as a sample of what your submission should look like. The file stage_2_detailed_class_info.csv contains detailed information about the positive and negative classes in the training set, and may be used to build more nuanced models.
- File descriptions
 - stage_2_train.csv the training set. Contains patientlds and bounding box / target information.
 - o **stage_2_sample_submission.csv** a sample submission file in the correct format. Contains patientlds for the test set. Note that the sample submission contains one box per image, but there is no limit to the number of bounding boxes that can be assigned to a given image.

- o **stage_2_detailed_class_info.csv** provides detailed information about the type of positive or negative class for each image.
- o **train_images** Contains DICOM images used to train the model.
- o test_images Contains DICOM images used to test the model.

Data fields

- o **patientId** A patientId. Each patientId corresponds to a unique image.
- o **x** the upper-left x coordinate of the bounding box.
- o **y** the upper-left y coordinate of the bounding box.
- o width the width of the bounding box.
- o **height** the height of the bounding box.
- o **Target** the binary Target, indicating whether this sample has evidence of pneumonia.

Implementation and EDA:

Importing Packages necessary for the project

import pandas as pd

import numpy as np

import os

import cv2

import seaborn as sns

import matplotlib.pyplot as plt

import tensorflow as tf

from tensorflow import keras

import PIL

from sklearn.model_selection import train_test_split

from sklearn.preprocessing import LabelEncoder

from sklearn.preprocessing import LabelBinarizer

from sklearn.preprocessing import StandardScaler

from tensorflow.keras.layers import Layer, Convolution2D, Flatten, Dense

from tensorflow.keras.layers import Concatenate, UpSampling2D, Conv2D, Reshape,

GlobalAveragePooling2D

from tensorflow.keras import losses,optimizers

from tensorflow.keras.layers import Dense, Activation,

Flatten, Dropout, MaxPooling 2D, Batch Normalization

from tensorflow.keras.models import Model from tensorflow.keras.models import Sequential from tensorflow.keras import backend as K import tensorflow.keras.utils as pltUtil from tensorflow.keras.utils import Sequence import math

import scipy.stats as stats

from sklearn.metrics import confusion_matrix

from matplotlib.patches import Rectangle

from tensorflow.keras import layers, models

from tensorflow.keras.utils import to_categorical

from tensorflow.keras.models import Sequential

from tensorflow.keras.layers import Flatten, Dense ,LeakyReLU,GlobalMaxPooling2D

from tensorflow.keras import regularizers, optimizers

from sklearn.metrics import r2_score

from tensorflow.keras.models import load_model

from tensorflow.keras.optimizers import RMSprop

 $from\ tensor flow. keras. preprocessing. image\ import\ Image Data Generator$

from tensorflow.keras.optimizers import Adam

from keras import applications

from tensorflow.keras.applications.vgg16 import VGG16

from tensorflow.keras.applications import ResNet50V2

from tensorflow.keras.applications import DenseNet121

from tensorflow.keras import datasets

from tensorflow.keras.layers import BatchNormalization

from tabulate import tabulate

from texttable import Texttable

import warnings

warnings.filterwarnings('ignore')

EDA on stage_2_train_labels File

train_labels= pd.read_csv('Dataset - New/stage_2_train_labels.csv')
print('First five rows of Training set:\n', train_labels.head())

First five rows of Training set:

```
patientId x y width height Target

0 0004cfab-14fd-4e49-80ba-63a80b6bddd6 NaN NaN NaN NaN 0

1 00313ee0-9eaa-42f4-b0ab-c148ed3241cd NaN NaN NaN NaN 0

2 00322d4d-1c29-4943-afc9-b6754be640eb NaN NaN NaN NaN 0

3 003d8fa0-6bf1-40ed-b54c-ac657f8495c5 NaN NaN NaN NaN 0

4 00436515-870c-4b36-a041-de91049b9ab4 264.0 152.0 213.0 379.0
```

print(train_labels.iloc[0])

patientId 0004cfab-14fd-4e49-80ba-63a80b6bddd6

x NaN
y NaN
width NaN
height NaN
Target 0

Name: 0, dtype: object

train_labels.shape

(30227, 6)

The data frame has 30227 rows and 6 Columns

train_labels.size

181362

train_labels.describe

```
<br/>
<bound method NDFrame.describe of
                                              patientId x y width height \
   0004cfab-14fd-4e49-80ba-63a80b6bddd6 NaN NaN NaN
                                                        NaN
   00313ee0-9eaa-42f4-b0ab-c148ed3241cd NaN NaN NaN
                                                        NaN
2
   00322d4d-1c29-4943-afc9-b6754be640eb NaN NaN NaN
                                                        NaN
3
   003d8fa0-6bf1-40ed-b54c-ac657f8495c5 NaN NaN NaN NaN
   00436515-870c-4b36-a041-de91049b9ab4 264.0 152.0 213.0 379.0
4
30222 c1ec14ff-f6d7-4b38-b0cb-fe07041cbdc8 185.0 298.0 228.0 379.0
30223 c1edf42b-5958-47ff-a1e7-4f23d99583ba NaN NaN NaN NaN
30224 c1f6b555-2eb1-4231-98f6-50a963976431 NaN NaN NaN NaN
30225 c1f7889a-9ea9-4acb-b64c-b737c929599a 570.0 393.0 261.0 345.0
30226 c1f7889a-9ea9-4acb-b64c-b737c929599a 233.0 424.0 201.0 356.0
```

Target

0 0

```
1
     0
2
      0
3
      0
4
      1
30222
        1
30223
        0
30224
      0
30225
30226
      1
```

[30227 rows x 6 columns]>

train_labels['patientId'].nunique()


There are 26684 Unique patientlds in the data frame, the above shows that each patientld has multiple rows in the dataframe

```
# Let us find how many patients out of 26684 has Pneumonia or not in our dataset # Creating Variables with and without Pneumonia to final the percentage distribution of data one = train_labels[train_labels.Target == 1].drop_duplicates('patientId').shape[0] zero = train_labels[train_labels.Target == 0].drop_duplicates('patientId').shape[0] total = train_labels.drop_duplicates('patientId').shape[0] print(f'No of Patients with Pneumonia: {one} - {round(one/total*100, 0)}% of the entire dataset') print(f'No of Patients without Pneumonia: {zero} - {round(zero/total*100, 0)}% of the entire dataset')
```

_ = train_labels.drop_duplicates('patientId').drop_duplicates('patientId')['Target'].value_counts().plot(kind = 'pie', autopct = '%.0f%%', labels = ['Negative', 'Positive'], figsize = (10, 6))

No of Patients with Pneumonia: 6012 - 23.0% of the entire dataset

No of Patients without Pneumonia: 20672 - 77.0% of the entire dataset



Checking nulls in bounding box columns: print('Number of nulls in bounding box columns: {}'.format(train_labels[['x', 'y', 'width', 'height']].isnull().sum().to_dict()))

Number of nulls in bounding box columns: {'x': 20672, 'y': 20672, 'width': 20672, 'height': 20672} Thus, we can see that the number of nulls in bounding box columns are equal to the number of 0's we have in the Target column.

EDA on stage_2_detailed_class_info File

class_labels = pd.read_csv('Dataset - New/stage_2_detailed_class_info.csv')
print('First five rows of Class label dataset are:\n', class_labels.head())

First five rows of Class label dataset are:

patientId class

0 0004cfab-14fd-4e49-80ba-63a80b6bddd6 No Lung Opacity / Not Normal

1 00313ee0-9eaa-42f4-b0ab-c148ed3241cd No Lung Opacity / Not Normal

2 00322d4d-1c29-4943-afc9-b6754be640eb No Lung Opacity / Not Normal

3 003d8fa0-6bf1-40ed-b54c-ac657f8495c5 Normal

4 00436515-870c-4b36-a041-de91049b9ab4 Lung Opacity

Some information about the data field present in the 'stage_2_detailed_class_info.csv' are: patientId - A patientId. Each patientId corresponds to a unique image class - Have three values depending on what is the current state of the patient's lung: 'No Lung Opacity / Not Normal', 'Normal' and 'Lung Opacity'.

Shape of the dataframe class_labels.shape

(30227, 2)

The data frame has 30227 rows and 2 Columns

Size of the dataframe class_labels.size 60454

Basic info about the file class_labels.describe

<bc< th=""><th>und method NDFrame.describe of</th><th>patientId</th><th>class</th></bc<>	und method NDFrame.describe of	patientId	class									
0	0004cfab-14fd-4e49-80ba-63a80b6bddd6	No Lung Opacity / Not Normal										
1	1 00313ee0-9eaa-42f4-b0ab-c148ed3241cd No Lung Opacity / Not Normal											
2	00322d4d-1c29-4943-afc9-b6754be640eb	No Lung Opacity / Not Normal										
3	003d8fa0-6bf1-40ed-b54c-ac657f8495c5	Normal										
4	00436515-870c-4b36-a041-de91049b9ab4	Lung Opacity										
• • •												
302	22 c1ec14ff-f6d7-4b38-b0cb-fe07041cbdc8	Lung Opacity										
302	23 c1edf42b-5958-47ff-a1e7-4f23d99583b	a Normal										
302	24 c1f6b555-2eb1-4231-98f6-50a96397643	Normal										
302	25 c1f7889a-9ea9-4acb-b64c-b737c92959	Pa Lung Opacity										
302	226 c1f7889a-9ea9-4acb-b64c-b737c92959	Pa Lung Opacity										
[30.	227 rows x 2 columns]>											

class_labels['patientId'].nunique()

26684

Printing the unique values in Class column to do the further analysis print(sorted(class_labels['class'].unique()))

['Lung Opacity', 'No Lung Opacity / Not Normal', 'Normal']

There are 26684 Unique patientlds in the data frame, the above shows that each patientld has multiple rows in the dataframe

The data is distributed across in the class column as below (class_labels['class'].value_counts() / class_labels['patientId'].count()) * 100


No Lung Opacity / Not Normal 39.107421

Lung Opacity 31.610812 Normal 29.281768

Name: class, dtype: float64

figsize = (10, 6)

_ = class_labels['class'].value_counts().sort_index(ascending = False).plot(kind = 'pie', autopct = '%.0f%%').set_ylabel(")



Checking nulls in class_labels:

print('Number of nulls in class columns: {}'.format(class_labels['class'].isnull().sum()))

Number of nulls in class columns: 0

Checking whether each patientId has only one type of class or not class_labels.groupby(['patientId'])['class'].nunique().max()

1

So, each patient is associated with 1 class

Merging of both the dataframes

Both train_lables and class_labels can be merged because of the same number of rows and the same number of unique patient IDs available in both the dataframes.

Merging the two dataset - 'train_labels' and 'class_labels':
training_data = pd.concat([train_labels, class_labels['class']], axis = 1)
print('After merging, the dataset looks like: \n')
training_data.head()

After merging, the dataset looks like:

	patientId	Χ	У	wi dth	hei ght	Tar get	class
0	0004cfab-14fd-4e49-80 ba-63a80b6bddd6	N a N	N a N	Na N	Na N	0	No Lung Opacity / Not Normal
1	00313ee0-9eaa-42f4-b0 ab-c148ed3241cd	N a N	N a N	Na N	Na N	0	No Lung Opacity / Not Normal
2	00322d4d-1c29-4943-af c9-b6754be640eb	N a N	N a N	Na N	Na N	0	No Lung Opacity / Not Normal
3	003d8fa0-6bf1-40ed-b5 4c-ac657f8495c5	N a N	N a N	Na N	Na N	0	Normal
4	00436515-870c-4b36-a 041-de91049b9ab4	26 4. 0	15 2. 0	21 3.0	379 .0	1	Lung Opacity

Shape of the dataframe training_data.shape

(30227, 7)

Size of the dataframe

training_data.size

211589

basic info

training_data.describe

```
<br/>
<bound method NDFrame.describe of
                                            patientId x y width height \
   0004cfab-14fd-4e49-80ba-63a80b6bddd6 NaN NaN NaN
                                                        NaN
   00313ee0-9eaa-42f4-b0ab-c148ed3241cd NaN NaN NaN
                                                        NaN
   00322d4d-1c29-4943-afc9-b6754be640eb NaN NaN NaN
2
                                                        NaN
3
   003d8fa0-6bf1-40ed-b54c-ac657f8495c5 NaN NaN NaN NaN
   00436515-870c-4b36-a041-de91049b9ab4 264.0 152.0 213.0 379.0
4
30222 c1ec14ff-f6d7-4b38-b0cb-fe07041cbdc8 185.0 298.0 228.0 379.0
30223 c1edf42b-5958-47ff-a1e7-4f23d99583ba NaN NaN NaN NaN
30224 c1f6b555-2eb1-4231-98f6-50a963976431 NaN NaN NaN NaN
30225 c1f7889a-9ea9-4acb-b64c-b737c929599a 570.0 393.0 261.0 345.0
30226 c1f7889a-9ea9-4acb-b64c-b737c929599a 233.0 424.0 201.0 356.0
```

Tai	rget	class					
0	0 No	Lung Opacity / Not Normal					
1	0 No	Lung Opacity / Not Normal					
2	0 No	Lung Opacity / Not Normal					
3	0	Normal					
4	1 Lung Opacity						
30222	1	Lung Opacity					
30223	0	Normal					
30224	0	Normal					
30225	1	Lung Opacity					
30226	1	Lung Opacity					

[30227 rows x 7 columns]>

training_data.isna().sum() #checking for null values

 patientId
 0

 x
 20672

 y
 20672

 width
 20672

 height
 20672

 Target
 0

 class
 0

 dtype: int64

training_data['Target'].value_counts() #Value count of the target

0 206721 9555


Name: Target, dtype: int64

training_data['class'].value_counts() #Frequency of the class


No Lung Opacity / Not Normal 11821

Lung Opacity 9555
Normal 8851
Name: class, dtype: int64

pd.value_counts(training_data['Target']).plot(kind="bar") #Bar plot



pd.value_counts(training_data['class']).plot(kind="bar") #Bar plot



training_data.head(5)

patien tId	X	у	wi dth	hei ght	Tar get	cl a s	
0	0004cfab-14fd-4e49-80ba -63a80b6bddd6	N a N	Na N	Na N	Na N	0	No Lung Opacity / Not Normal
1	00313ee0-9eaa-42f4-b0a b-c148ed3241cd	N a N	Na N	Na N	Na N	0	No Lung Opacity / Not Normal
2	00322d4d-1c29-4943-afc 9-b6754be640eb	N a N	Na N	Na N	Na N	0	No Lung Opacity / Not Normal


Let us now see how the target in train_labels and class in class_labels are associated fig, ax = plt.subplots(nrows = 1, figsize = (12, 6))

temp = training_data.groupby('Target')['class'].value_counts()

data_target_class = pd.DataFrame(data = {'Values': temp.values}, index =
temp.index).reset_index()

sns.barplot(ax = ax, x = 'Target', y = 'Values', hue = 'class', data = data_target_class, palette = 'Set3')

plt.title('Class and Target for Chest Exams')




So, from the above we can conclude the below Grouping:

If the Target is 0, the class would be either - No Lung Opacity / Not Normal or Normal

If the Target is 1, the class would be Lung Opacity

#Distribution Plot
target1 = training_data[training_data['Target']==1]
sns.set_style('whitegrid')
plt.figure()
fig, ax = plt.subplots(2,2,figsize=(10,10))

sns.distplot(target1['x'],kde=True,bins=50, color="red", ax=ax[0,0]) sns.distplot(target1['y'],kde=True,bins=50, color="blue", ax=ax[0,1]) sns.distplot(target1['width'],kde=True,bins=50, color="green", ax=ax[1,0]) sns.distplot(target1['height'],kde=True,bins=50, color="yellow", ax=ax[1,1]) locs, labels = plt.xticks() plt.tick_params(axis='both', which='major', labelsize=10) plt.show()



#Heatmap

fig, ax = plt.subplots(1,1,figsize=(10,10))

target2 = target1.sample(1000)

target2['xc'] = target2['x'] + target2['width'] / 2


target2['yc'] = target2['y'] + target2['height'] / 2

plt.title("Points of Lung Opacity: red rectangles over yellow rectangles:")

target2.plot.scatter(x='xc', y='yc', xlim=(0,1024), ylim=(0,1024), ax=ax, alpha=0.8, marker=".", color="red")

for i, crt_sample in target2.iterrows():

ax.add_patch(Rectangle(xy=(crt_sample['x'],
crt_sample['y']),width=crt_sample['width'],height=crt_sample['height'],alpha=3.5e-3,
color="yellow"))
plt.show()



Train Folder Images
print('Number of images in training images folders are: {}.'.format(len(os.listdir('Dataset New/stage_2_train_images (1)-001'))))

Number of images in training images folders are: 26684.

We can see that in the training images folder we have just 26684 images which is the same as that of the unique patientld's present in either of the csv files. Thus, we can say that each of the unique patientld's present in either of the csv files corresponds to an image present in the folder.

Creating images dataframe with 2 columns - path of the dicom file and associated patientId # Every DCM filename is the patientId from glob import glob training_image_path = 'Dataset - New/stage_2_train_images (1)-001' images = pd.DataFrame({'path': glob(os.path.join(training_image_path, '*.dcm'))}) images['patientId'] = images['path'].map(lambda x:os.path.splitext(os.path.basename(x))[0])

print('Columns in the training images dataframe: {}'.format(list(images.columns)))

Columns in the training images dataframe: ['path', 'patientId']

Merging the images data frame with training_data dataframe training_data = training_data.merge(images, on = 'patientId', how = 'left') print('After merging the two dataframe, the training_data has {} rows and {} columns.'.format(training_data.shape[0], training_data.shape[1]))

After merging the two dataframe, the training_data has 30227 rows and 8 columns.

training_data.head()

	patientId	Χ	у	wi dth	hei ght	Tar get	class	path
0	0004cfab-14f d-4e49-80ba-6 3a80b6bddd6	N a N	N a N	Na N	Na N	0	No Lung Opacit y / Not Normal	Dataset - New/stage_2_train_ima ges (1)-001\000
1	00313ee0-9ea a-42f4-b0ab-c 148ed3241cd	N a N	N a N	Na N	Na N	0	No Lung Opacit y / Not Normal	Dataset - New/stage_2_train_ima ges (1)-001\003
2	00322d4d-1c 29-4943-afc9- b6754be640e b	N a N	N a N	Na N	Na N	0	No Lung Opacit y / Not Normal	Dataset - New/stage_2_train_ima ges (1)-001\003
3	003d8fa0-6bf 1-40ed-b54c-a c657f8495c5	N a N	N a N	Na N	Na N	0	Normal	Dataset - New/stage_2_train_ima ges (1)-001\003
4	00436515-87 0c-4b36-a041- de91049b9ab 4	26 4. 0	15 2. 0	21 3.0	379 .0	1	Lung Opacit y	Dataset - New/stage_2_train_ima ges (1)-001\004

Shape of the dataframe training_data.shape

(30227, 8)

Size of the dataframe training_data.size

241816

basic info of the dataframe training_data.describe

```
<bound method NDFrame.describe of</p>
                                              patientId
                                                         x y width height \
   0004cfab-14fd-4e49-80ba-63a80b6bddd6 NaN NaN NaN
0
                                                         NaN
   00313ee0-9eaa-42f4-b0ab-c148ed3241cd NaN
1
                                             NaN
                                                   NaN
                                                         NaN
   00322d4d-1c29-4943-afc9-b6754be640eb NaN NaN NaN
                                                         NaN
3
   003d8fa0-6bf1-40ed-b54c-ac657f8495c5 NaN NaN NaN
                                                        NaN
4
   00436515-870c-4b36-a041-de91049b9ab4 264.0 152.0 213.0 379.0
30222 c1ec14ff-f6d7-4b38-b0cb-fe07041cbdc8 185.0 298.0 228.0 379.0
30223 c1edf42b-5958-47ff-a1e7-4f23d99583ba NaN NaN NaN
                                                         NaN
30224 c1f6b555-2eb1-4231-98f6-50a963976431 NaN NaN NaN NaN
30225 c1f7889a-9ea9-4acb-b64c-b737c929599a 570.0 393.0 261.0 345.0
30226 c1f7889a-9ea9-4acb-b64c-b737c929599a 233.0 424.0 201.0 356.0
```

```
class \
   Target
      0 No Lung Opacity / Not Normal
0
1
      0 No Lung Opacity / Not Normal
2
      0 No Lung Opacity / Not Normal
3
                    Normal
      1
                 Lung Opacity
30222
                   Lung Opacity
                      Normal
30223
         0
30224
         0
                      Normal
30225
         1
                   Lung Opacity
30226
         1
                   Lung Opacity
```

```
path

Dataset - New/stage_2_train_images (1)-001\000...

Dataset - New/stage_2_train_images (1)-001\003...

Dataset - New/stage_2_train_images (1)-001\003...

Dataset - New/stage_2_train_images (1)-001\003...

Dataset - New/stage_2_train_images (1)-001\004...

...

30222 Dataset - New/stage_2_train_images (1)-001\c1e...
30223 Dataset - New/stage_2_train_images (1)-001\c1e...
30224 Dataset - New/stage_2_train_images (1)-001\c1f...
30225 Dataset - New/stage_2_train_images (1)-001\c1f...
```

30226 Dataset - New/stage_2_train_images (1)-001\c1f...

[30227 rows x 8 columns]>

Function to show images from train folder with bounding boxes

from matplotlib.patches import Rectangle import pydicom as dcm

```
def show_dicom_images(data, df, img_path):
 img_data = list(data.T.to_dict().values())
f, ax = plt.subplots(3, 3, figsize = (16, 18))
 for i, row in enumerate(img_data):
  image = row['patientId'] + '.dcm'
  path = os.path.join(img_path, image)
  data = dcm.read_file(path)
  rows = df[df['patientId'] == row['patientId']]
  data_img = dcm.dcmread(path)
  ax[i//3, i%3].imshow(data_img.pixel_array, cmap = plt.cm.bone)
  ax[i//3, i%3].axis('off')
  ax[i//3, i%3].set_title('ID: {}\nTarget: {}, Class: {}'\
                .format(row['patientId'], row['Target'],
                     row['class'], row['x'],
                     row['y'], row['width'],
                     row['height']))
  box_data = list(rows.T.to_dict().values())
  for j, row in enumerate(box_data):
```

Showing images without pneumonia

show_dicom_images(data = training_data.loc[(training_data['Target'] == 0)].sample(9), df = training_data, img_path = 'Dataset - New/stage_2_train_images (1)-001')

ID: 484eb1ad-b27e-4ff7-84e4-402753501616 Target: 0, Class: No Lung Opacity / Not Normal



ID: ea616acd-1b00-4342-94c5-7ad631ae9016 Target: 0, Class: No Lung Opacity / Not Normal



ID: eb832732-4c11-4126-bb10-f95eb61608bc Target: 0, Class: Normal



ID: d9eb090f-c64f-48a1-b38d-519d4dad2831 Target: 0, Class: Normal



ID: 57c51bb2-7795-48a1-8f14-e16e7c0ee3eb Target: 0, Class: Normal



ID: 91de1105-565e-4f72-8c39-7bb5a1246398 Target: 0, Class: Normal



ID: 8956c7cd-6593-48c3-a140-ba1be9fff0b8 Target: 0, Class: No Lung Opacity / Not Normal



ID: e4893d4a-f3d7-4d25-b733-da36f13f9981 Target: 0, Class: No Lung Opacity / Not Normal



ID: c848c8e3-50ac-4185-8317-e8baecaabc0e Target: 0, Class: Normal



Showing images with pneumonia

show_dicom_images(data = training_data.loc[(training_data['Target'] == 1)].sample(9), df = training_data, img_path = 'Dataset - New/stage_2_train_images (1)-001')

ID: 23815776-6d93-4835-867a-1b251b10b1d3 Target: 1, Class: Lung Opacity




ID: 826361cd-08ae-49c4-b535-25a4f14026f8 Target: 1, Class: Lung Opacity



ID: 8cddc2e8-0283-4be1-8fc4-b683cef00bbe Target: 1, Class: Lung Opacity



ID: 7bf09d53-83de-418a-ac0f-1c29a9dcf6fe Target: 1, Class: Lung Opacity




ID: becced78-9df7-4e9a-94c8-7cb750d41507 Target: 1, Class: Lung Opacity




ID: b16bcbb6-570b-450b-a9e5-b1e3663205bb Target: 1, Class: Lung Opacity



ID: 51b0e89f-608a-4040-9d0d-694ebe072c28 Target: 1, Class: Lung Opacity



ID: ba87678d-df8e-42eb-9a2b-3f8a6b34722a Target: 1, Class: Lung Opacity



ID: c5ba46fb-eddc-4824-b779-0620b0cd855c Target: 1, Class: Lung Opacity



EDA on Metadata of DICOM Images

Medical images are stored in a special format known as DICOM files (*.dcm). They contain a combination of header metadata as well as underlying raw image arrays for pixel data. In Python, one popular library to access and manipulate DICOM files is the pydicom module. To use the pydicom library, first find the DICOM file for a given patientld by simply looking for the matching file in the stage_2_train_images/ folder, and the use the pydicom.read_file() method to load the data

sample_patientId = train_labels['patientId'][0]
dcm_file = training_image_path+'/{}.dcm'.format(sample_patientId)
dcm_data = dcm.read_file(dcm_file)
print('Metadata of the image consists of \n', dcm_data)

```
Metadata of the image consists of
Dataset.file_meta
(0002, 0000) File Meta Information Group Length UL: 202
(0002, 0001) File Meta Information Version
                                            OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID
                                             UI: Secondary Capture Image Storage
(0002, 0003) Media Storage SOP Instance UID
1.2.276.0.7230010.3.1.4.8323329.28530.1517874485.775526
                                        UI: JPEG Baseline (Process 1)
(0002, 0010) Transfer Syntax UID
(0002, 0012) Implementation Class UID
                                           UI: 1.2.276.0.7230010.3.0.3.6.0
(0002, 0013) Implementation Version Name
                                              SH: 'OFFIS_DCMTK_360'
(0008, 0005) Specific Character Set
                                         CS: 'ISO_IR 100'
                                       UI: Secondary Capture Image Storage
(0008, 0016) SOP Class UID
(0008, 0018) SOP Instance UID
1.2.276.0.7230010.3.1.4.8323329.28530.1517874485.775526
(0008, 0020) Study Date
                                     DA: '19010101'
(0008, 0030) Study Time
                                     TM: '000000.00'
                                         SH: "
(0008, 0050) Accession Number
(0008, 0060) Modality
                                    CS: 'CR'
(0008, 0064) Conversion Type
                                       CS: 'WSD'
(0008, 0090) Referring Physician's Name
                                            PN: "
(0008, 103e) Series Description
                                       LO: 'view: PA'
(0010, 0010) Patient's Name
                                      PN: '0004cfab-14fd-4e49-80ba-63a80b6bddd6'
(0010, 0020) Patient ID
                                    LO: '0004cfab-14fd-4e49-80ba-63a80b6bddd6'
(0010, 0030) Patient's Birth Date
                                       DA: "
(0010, 0040) Patient's Sex
                                     CS: 'F'
(0010, 1010) Patient's Age
                                     AS: '51'
```

```
(0018, 0015) Body Part Examined
                                       CS: 'CHEST'
(0018, 5101) View Position
                                    CS: 'PA'
(0020, 000d) Study Instance UID
                                      UI:
1.2.276.0.7230010.3.1.2.8323329.28530.1517874485.775525
(0020, 000e) Series Instance UID
                                      UI:
1.2.276.0.7230010.3.1.3.8323329.28530.1517874485.775524
                                  SH: "
(0020, 0010) Study ID
(0020, 0011) Series Number
                                     IS: '1'
(0020, 0013) Instance Number
                                      IS: '1'
                                     CS: "
(0020, 0020) Patient Orientation
(0028, 0002) Samples per Pixel
                                      US: 1
(0028, 0004) Photometric Interpretation CS: 'MONOCHROME2'
(0028, 0010) Rows
                                 US: 1024
(0028, 0011) Columns
                                   US: 1024
(0028, 0030) Pixel Spacing
                                   DS: [0.14300000000000002, 0.143000000000000002]
(0028, 0100) Bits Allocated
                                   US: 8
(0028, 0101) Bits Stored
                                  US: 8
(0028, 0102) High Bit
                                 US: 7
(0028, 0103) Pixel Representation
                                      US: 0
                                        CS: '01'
(0028, 2110) Lossy Image Compression
(7fe0, 0010) Pixel Data
                                  OB: Array of 142006 elements
## Adding necessary Columns for EDA from Metadata of the DICOM Images
columns_to_add = ['Modality', 'PatientAge', 'PatientSex', 'BodyPartExamined', 'ViewPosition',
'ConversionType', 'Rows', 'Columns', 'PixelSpacing']
from tqdm import tqdm, tqdm_notebook
def parse_dicom_data(data_df, data_path):
for col in columns_to_add:
  data_df[col] = None
 image_names = os.listdir('Dataset - New/stage_2_train_images (1)-001/')
 for i, img_name in tqdm_notebook(enumerate(image_names)):
  imagepath = os.path.join('Dataset - New/stage_2_train_images (1)-001/', img_name)
  data_img = dcm.read_file(imagepath)
  idx = (data_df['patientId'] == data_img.PatientID)
  data_df.loc[idx, 'Modality'] = data_img.Modality
  data_df.loc[idx, 'PatientAge'] = pd.to_numeric(data_img.PatientAge)
```

```
data_df.loc[idx, 'PatientSex'] = data_img.PatientSex
data_df.loc[idx, 'BodyPartExamined'] = data_img.BodyPartExamined
data_df.loc[idx, 'ViewPosition'] = data_img.ViewPosition
data_df.loc[idx, 'ConversionType'] = data_img.ConversionType
data_df.loc[idx, 'Rows'] = data_img.Rows
data_df.loc[idx, 'Columns'] = data_img.Columns
data_df.loc[idx, 'PixelSpacing'] = str.format("{:4.3f}", data_img.PixelSpacing[0])
```

Parse DICOM Images to pull the metadata from them parse_dicom_data(training_data,'Dataset - New/stage_2_train_images (1)-001/')

print('So after parsing the information from the dicom images, our training_data data frame has {} rows and {} columns and it looks like:\n'.format(training_data.shape[0], training_data.shape[1])) training_data.head()

patientl	d :	k y	width	height	Target	class	path	Modality	PatientAge	PatientSex	BodyPartExamined	ViewPosition	${\sf Conversion Type}$	Rows	Columns	PixelSpacing
0004cfab 14fd-4e49 80ba 63a80b6bddd)- Nai	N NaN	NaN	NaN	0	No Lung Opacity / Not Normal	Dataset - New/stage_2_train_images (1)-001\000	CR	51	F	CHEST	РД	WSD	1024	1024	0.143
00313ee0 9eaa-42f4 b0ab c148ed3241c	l- Nai	N NaN	NaN	NaN	0	No Lung Opacity / Not Normal	Dataset - New/stage_2_train_images (1)-001\003	CR	48	F	CHEST	PA	WSD	1024	1024	0.194
00322d4c 1c29-4943 afc9 b6754be640e	}- - Naf	N NaN	NaN	NaN	0	No Lung Opacity / Not Normal	Dataset - New/stage_2_train_images (1)-001\003	CR	19	М	CHEST	ДР	WSD	1024	1024	0.168
003d8fa0 6bf1-40ec b54 ac657f8495d	l- Nai	N NaN	NaN	NaN	0	Normal	Dataset - New/stage_2_train_images (1)-001\003	CR	28	М	CHEST	PA	WSD	1024	1024	0.143
00436511 870c-4b36 a04 de91049b9ab	5- 264.	0 152.0	213.0	379.0	1	Lung Opacity	Dataset - New/stage_2_train_images (1)-001\004	CR	32	F	CHEST	АР	WSD	1024	1024	0.139

EDA on Modality

print('Modality for the images obtained is: {} \n'.format(training_data['Modality'].unique()[0])) Modality for the images obtained is: CR

EDA on Body Part Examined

print('The images obtained are of {} areas.'.format(training_data['BodyPartExamined'].unique()[0]))
The images obtained are of CHEST areas.

EDA on Understanding Different Positions

(training_data['ViewPosition'].value_counts() / training_data['patientId'].count()) * 100

AP 50.607073 PA 49.392927 Name: ViewPosition, dtype: float64

AP : Anterior/Posterior (AP)
PA : Posterior/Anterior (PA)

EDA on Conversion Type

print('Conversion Type for the data in Training Data: ',
training_data['ConversionType'].unique()[0])

Conversion Type for the data in Training Data: WSD

EDA on Gender

```
(training_data['PatientSex'].value_counts() / training_data['patientId'].count()) * 100

M 56.955702

F 43.044298

Name: PatientSex, dtype: float64
```


Function to plot graphs on different attributes of the dataframe

```
def drawgraphs(data_file, columns, hue = False, width = 15, showdistribution = True):
 if (hue):
  print('Creating graph for: {} and {}'.format(columns, hue))
  print('Creating graph for : {}'.format(columns))
 length = len(columns) * 6
total = float(len(data_file))
fig, axes = plt.subplots(nrows = len(columns) if len(columns) > 1 else 1, ncols = 1, figsize =
(width, length))
 for index, content in enumerate(columns):
  plt.title(content)
  currentaxes = 0
  if (len(columns) > 1):
   currentaxes = axes[index]
  else:
   currentaxes = axes
  if (hue):
   sns.countplot(x = columns[index], data = data_file, ax = currentaxes, hue = hue)
  else:
   sns.countplot(x = columns[index], data = data_file, ax = currentaxes)
  if(showdistribution):
   for p in (currentaxes.patches):
    height = p.get_height()
```

if (height > 0 and total > 0): currentaxes.text(p.get_x() + p.get_width()/2., height + 3, $\frac{1.2f}{1.2f}$ format(100*height/total), ha = "center")


EDA on PatientSex and Count of patientId - Hue by Class

drawgraphs(data_file = training_data.drop_duplicates('patientId'), columns = ['PatientSex'], hue = 'class', width = 20, showdistribution = True)



EDA on PatientSex and Count of patientId. Split by Target

drawgraphs(data_file = training_data.drop_duplicates('patientId'), columns = ['PatientSex'], hue = 'Target', width = 20, showdistribution = True)



EDA on PatientAge


Split PatientAge based on into PatientAgeBins custom_array = np.linspace(0, 100, 11) training_data['PatientAgeBins'] = pd.cut(training_data['PatientAge'], custom_array) training_data.drop_duplicates('patientId')['PatientAgeBins'].value_counts()

(50.0, 60.0] 6619 (40.0, 50.0] 5048 (60.0, 70.0] 4279 (30.0, 40.0]3936 (20.0, 30.0]3205 (10.0, 20.0]1468 (70.0, 80.0] 1452 (0.0, 10.0]435 (80.0, 90.0] 233 (90.0, 100.0]

Name: PatientAgeBins, dtype: int64


Distribution on Count of Patients based on PatientAgeBins

drawgraphs(data_file = training_data.drop_duplicates('patientId'), columns = ['PatientAgeBins'], width = 20, showdistribution = True)




As per our dataset the major images belong to age group between 50 and 60

drawgraphs(data_file = training_data.drop_duplicates('patientId'), columns = ['PatientAgeBins'], hue = 'PatientSex', width = 20, showdistribution = True)




The Images for Male are more than of Female

drawgraphs(data_file = training_data.drop_duplicates('patientId'), columns = ['PatientAgeBins'], hue = 'class', width = 20, showdistribution = True)



drawgraphs(data_file = training_data.drop_duplicates('patientId'), columns = ['PatientAgeBins'], hue = 'Target', width = 20, showdistribution = True)



Test Folder Images

print('Number of images in testing images folders are: {}.'.format(len(os.listdir('Dataset - New/stage_2_test_images'))))

Number of images in testing images folders are: 3000.

Image Classification

Image classification is the process of predicting a specific class, or label, for something that is defined by a set of data points. In this case, we are building a deep learning architecture to classify the given patient's chest X-ray image as a case of pneumonia or not. The dataset primarily was not designed for image classification and this is where we were creative. Usually, the data is organized as image files in folders for each class. The generator will take each

folder as the classes to be trained. Keras generator function flow_from_dataframe() allows us to organize the data as a dataframe where we have the image filename and their corresponding target. Here we only use the train images as it is labeled. The DICOM images are converted into JPEG images as it is easier for the model to process and gives better metrics. Then, pre-processing of the converted images is done after splitting the data into train, validation and test sets. The Keras ImageDataGenerator takes in these images, transforms them with the arguments we give and feeds them into the model, which results in a better model performance, although it takes more time.

Transfer learning allows us to use a pre-existing model, trained on a huge dataset, for our own tasks. Consequently reducing the cost of training new deep learning models and since the datasets have been vetted, we can be assured of the quality. We have built three image classification models with transfer learning:

- 1.VGG16 Model
- 2.ResNet50 Model
- 3.DenseNet121 Model

The three models are fairly consistent and have a similar framework. In all three models, pre-trained imagenet weights are used as it provides a standard measure for model comparison and helps in model coverage with less number of epochs. The input shape is also the same and include_top =False as we are using pre-trained material and we do not want the final dense layer that gives classification prediction. After the architecture, we compile the model with three metrics: accuracy, precision and recall with loss being set to categorical cross entropy and binary cross entropy as it is a classification problem and the optimizer is Adam. We then display the summary. Then, it is fit with the training and validation data along with other parameters. The average accuracy of the training data is shown and plots displaying the accuracy and loss metrics against the number epochs are visualized. Finally, the test metrics are displayed individually and in a summary table.

VGG-16 Model:

We set some layers frozen. To "freeze" a layer means to exclude it from training. Flattening of the last layer is done. It converts a multi-dimensional feature map to a single dimension. The prediction Dense layer has 2 nodes and softmax activation.

ResNet50 Model and Dense121 Model:

Even in the ResNet50,we freeze some layers. The ResNet50 model and DenseNet121 model architectures are almost identical. Dropout threshold is set to 0.5 to prevent overfitting and to regularize the model. There are 2 neurons in the output layer with sigmoid activation.

In this context, accuracy is the most important metric for image classification, not taking anything away from the other metrics. The train accuracy is most important as we are interested in the performance of the model in the real world setting. All three models have accuracies in the range of 70% - 80% with VGG16 having the best test accuracy. The ResNet50 model closely follows and DenseNet121 has the least test accuracy. We can add more trainable layers for performance improvements.

Converting the train images from DICOM to JPEG

```
import pydicom as dicom
PNG = False
jpg_folder_path = "Dataset - New/JPG_test" #Output folder
images_path = os.listdir(training_image_path)
for n, image in enumerate(images_path):
  ds = dicom.dcmread(os.path.join(training_image_path, image))
  pixel_array_numpy = ds.pixel_array
  if PNG == False:
    image = image.replace('.dcm', '.jpg')
  else:
    image = image.replace('.dcm', '.png')
  cv2.imwrite(os.path.join(jpg_folder_path, image), pixel_array_numpy)
  if n \% 50 == 0:
    print('{} image converted'.format(n))
0 image converted
50 image converted
100 image converted
150 image converted
200 image converted
26500 image converted
26550 image converted
26600 image converted
26650 image converted
```

```
df = pd.merge(train_labels, class_labels, on="patientId")
#Merging the images with their labels in a dataframe
images = pd.DataFrame({'path': glob(os.path.join(jpg_folder_path, '*.jpg'))})
images['patientId'] = images['path'].map(lambda x:os.path.splitext(os.path.basename(x))[0])
print('Columns in the training images dataframe: {}'.format(list(images.columns)))
merge = df.merge(images, on = 'patientId', how = 'left')
print('After merging the two dataframe, the training_data has {} rows and {}
columns.'.format(df.shape[0], df.shape[1]))
Columns in the training images dataframe: ['path', 'patientId']
After merging the two dataframe, the training_data has 37629 rows and 7 columns.
#Changing the target variable to string as it a classification problem
new['Target'] = new['Target'].astype(str)
#Splitting the dataset into train, validation and test with the below size
t1=new.iloc[0:31625] #train
t2=new.iloc[31625:32625] #validation
t3=new.iloc[32625:37625] #test
print("t1 : ", t1.shape )
print("t2:", t2.shape)
print("t3:", t3.shape)
t1: (31625, 2)
t2: (1000, 2)
t3: (5000, 2)
```

Image Pre-processing


return score

```
from keras.preprocessing.image import ImageDataGenerator
train_datagenerator = ImageDataGenerator(rescale=1. /
255,shear_range=0.2,zoom_range=0.2,rotation_range=20,brightness_range=(1.2,
1.5),horizontal_flip=True)
test_datagenerator = ImageDataGenerator(rescale=1. / 255)
val_datagenerator = ImageDataGenerator(rescale=1. / 255)
train_generator =
train_datagenerator.flow_from_dataframe(dataframe=t1,x_col='path',y_col="Target",target_size
=(224, 224),batch_size=16,class_mode='categorical')
val_generator =
val_datagenerator.flow_from_dataframe(dataframe=t2,x_col='path',y_col="Target",target_size=(
224, 224),batch_size=16,class_mode='categorical')
test_generator =
test_datagenerator.flow_from_dataframe(dataframe=t3,x_col='path',y_col="Target",target_size=
(224, 224),batch_size=16,class_mode='categorical')
Found 31625 validated image filenames belonging to 2 classes.
Found 1000 validated image filenames belonging to 2 classes.
Found 5000 validated image filenames belonging to 2 classes.
Functions for Re-usability:
headAccTable = ["Model","Train Accuracy", "Test Accuracy"]
headDataMetrics = ["Loss","Accuracy", "Precision", "Recall"]
def printTable(data, head):
   # display table
  print(tabulate(data, headers=head, tablefmt="fancy_grid"))
def getScore(modelHistory):
  score = np.mean(modelHistory.history['accuracy'])
  print('The Average Training Accuracy is', score)
```

```
def plotAccAndLoss(modelHistory):
    #Plots-accuracy and loss in train and validation data
    fig, ax = plt.subplots(1, 2, figsize=(10, 3))
    for i, met in enumerate(['accuracy', 'loss']):
        ax[i].plot(modelHistory.history[met])
        ax[i].plot(modelHistory.history['val_' + met])
        ax[i].set_title('Model {}'.format(met))
        ax[i].set_xlabel('epochs')
        ax[i].set_ylabel(met)
        ax[i].legend(['train', 'val'])
```

VGG-16 Model

A convolutional neural network is also known as a ConvNet, which is a kind of artificial neural network. A convolutional neural network has an input layer, an output layer, and various hidden layers. VGG16 is a type of CNN (Convolutional Neural Network) that is considered to be one of the best computer vision models to date. The creators of this model evaluated the networks and increased the depth using an architecture with very small (3×3) convolution filters, which showed a significant improvement on the prior-art configurations. They pushed the depth to 16-19 weight layers making it approx -138 trainable parameters. VGG16 is an object detection and classification algorithm which is able to classify 1000 images of 1000 different categories with 92.7% accuracy.



Vgg model building

```
vgg = VGG16(weights='imagenet',include_top=False,input_shape=(224,224,3))
for layer in vgg.layers:
    layer.trainable = False #layers as non-trainable
#flattening last layer
x = Flatten()(vgg.output)
#Dense layer
```

predictions = Dense(2,activation='softmax')(x)

model1 = Model(inputs=vgg.input, outputs=predictions)

model1.compile(loss='categorical_crossentropy',optimizer='adam',metrics=['accuracy','Precisio n','Recall'])

model1.summary()

Model: "model"

Layer (type) Output Shape Param # ______ [(None, 224, 224, 3)] 0 input_1 (InputLayer) block1_conv1 (Conv2D) (None, 224, 224, 64) 1792 block1_conv2 (Conv2D) (None, 224, 224, 64) 36928 block1_pool (MaxPooling2D) (None, 112, 112, 64) 0 block2_conv1 (Conv2D) (None, 112, 112, 128) 73856 block2_conv2 (Conv2D) (None, 112, 112, 128) 147584 block2_pool (MaxPooling2D) (None, 56, 56, 128) 0 block3_conv1 (Conv2D) (None, 56, 56, 256) 295168 block3_conv2 (Conv2D) (None, 56, 56, 256) 590080 block3_conv3 (Conv2D) (None, 56, 56, 256) 590080 block3_pool (MaxPooling2D) (None, 28, 28, 256) 0 block4_conv1 (Conv2D) (None, 28, 28, 512) 1180160 block4_conv2 (Conv2D) (None, 28, 28, 512) 2359808 block4_conv3 (Conv2D) (None, 28, 28, 512) 2359808 block4_pool (MaxPooling2D) (None, 14, 14, 512) block5_conv1 (Conv2D) (None, 14, 14, 512) 2359808 block5_conv2 (Conv2D) (None, 14, 14, 512) 2359808 block5_conv3 (Conv2D) (None, 14, 14, 512) 2359808 block5_pool (MaxPooling2D) (None, 7, 7, 512) (None, 25088) flatten (Flatten) 0 dense (Dense) (None, 2) 50178 ______

Total params: 14,764,866 Trainable params: 50,178

Non-trainable params: 14,714,688

history1 = model1.fit(trai

Epoch 1/10

$model 1. fit (train_generator, validation_data=val_generator, steps_per_epoch=25, epochs=10)$

Epoch 2/10

25/25 [==========================] - 226s 9s/step - loss: 0.8257 - accuracy: 0.6675 - precision: 0.6675 - recall: 0.6675 - val_loss: 0.6705 - val_accuracy: 0.6890 - val_precision: 0.6890 - val_recall: 0.6890

Epoch 3/10

25/25 [=====================] - 237s 10s/step - loss: 0.6190 - accuracy: 0.7455 - precision: 0.7455 - recall: 0.7455 - val_loss: 0.4521 - val_accuracy: 0.7880 - val_precision: 0.7880 - val_recall: 0.7880

Epoch 4/10

Epoch 5/10

Epoch 6/10

25/25 [======================] - 234s 10s/step - loss: 0.6016 - accuracy: 0.7300 - precision: 0.7300 - recall: 0.7300 - val_loss: 0.6002 - val_accuracy: 0.7310 - val_precision: 0.7310 - val_recall: 0.7310

Epoch 7/10

Epoch 8/10

Epoch 9/10


25/25 [==========================] - 250s 10s/step - loss: 0.6067 - accuracy: 0.7475 - precision: 0.7475 - recall: 0.7475 - val_loss: 0.9398 - val_accuracy: 0.6270 - val_precision: 0.6270 - val_recall: 0.6270

Epoch 10/10

#Training Accuracy of the VGG-16 Model score1 = getScore(history1)

The Average Training Accuracy is 0.7188047051429749

#Plots-accuracy and loss in train vs validation data plotAccAndLoss(history1)



result1 = model1.evaluate_generator(test_generator, 5000) printTable([result1], headDataMetrics)

Loss Accuracy Precision	Recall	
0.587347 0.754 0.754	0.754	

ResNet50 Model

ResNet stands for Residual Network. This model was immensely successful, as can be ascertained from the fact that its ensemble won the top position at the ILSVRC 2015 classification competition with an error of only 3.57%. ResNet has many variants that run on the same concept but have different numbers of layers. Resnet50 is used to denote the variant that can work with 50 neural network layers.

ResNet50 Model Architecture Zero Padding Input **Sonv Block** Conv Block **Batch Norm** Conv Block Conv Block Output Flattening Block Block Block Block **Avg Pool** CONV ReLu ₽ C ₽ Stage 1 Stage 2 Stage 3 Stage 4 Stage 5

Resnet50 Model building

```
resnet50 = ResNet50V2(weights = "imagenet", input_shape = (224,224,3), include_top = False) for layer in resnet50.layers:
```

layer.trainable = False

model2 = Sequential()

model2.add(resnet50)

model2.add(Flatten())

model2.add(Dense(units = 128, activation = "relu"))

model2.add(Dropout(0.5))

model2.add(Dense(units = 2, activation = "sigmoid"))

model2.compile(optimizer = 'adam', loss = 'binary_crossentropy', metrics = ['accuracy','Precision','Recall'])

model2.summary()

Model: "sequential"

Layer (type) Output Shape Param # ______ resnet50v2 (Functional) (None, 7, 7, 2048) 23564800 flatten_1 (Flatten) (None, 100352) 0 dense_1 (Dense) (None, 128) 12845184 dropout (Dropout) (None, 128) 0 dense_2 (Dense) (None, 2) 258 ______

Total params: 36,410,242 Trainable params: 12,845,442

```
history2 =
```

model2.fit(train_generator,validation_data=val_generator,steps_per_epoch=25,epochs=10)

Epoch 1/10

Epoch 2/10

Epoch 3/10

Epoch 4/10

Epoch 5/10

25/25 [====================] - 85s 3s/step - loss: 0.6649 - accuracy: 0.6600 - precision: 0.5994 - recall: 0.9275 - val_loss: 0.6127 - val_accuracy: 0.7660 - val_precision: 0.6063 - val_recall: 0.8240

Epoch 6/10

Epoch 7/10

25/25 [=====================] - 87s 4s/step - loss: 0.6241 - accuracy: 0.6350 - precision: 0.5919 - recall: 0.9100 - val_loss: 0.6726 - val_accuracy: 0.4450 - val_precision: 0.5820 - val_recall: 0.8130


Epoch 8/10

Epoch 9/10> ====] - 84s 3s/step - loss: 0.6497 - accuracy: 0.6725 - precision: 0.5890 - recall: 0.9100 - val_loss: 0.5576 - val_accuracy: 0.7630 - val_precision: 0.6850 - val_recall: 0.8090

#Training Accuracy of the ResNet50 Model score2 = getScore(history2)

The Average Training Accuracy is 0.6725000083446503

#Plots-accuracy and loss in train vs validation data plotAccAndLoss(history2)



#Test Data Metrics-Accuracy,Precision and Recall
result2 = model2.evaluate_generator(test_generator, 5000)
printTable([result2], headDataMetrics)


Loss Accuracy Precision Recall			
0.50237 0.8024 0.74935 0.8652			

DenseNet121 Model

DenseNet (Dense Convolutional Network) is an architecture that focuses on making the deep learning networks go even deeper, but at the same time making them more efficient to train, by using shorter connections between the layers. DenseNet is a convolutional neural network where each layer is connected to all other layers that are deeper in the network, that is, the first layer is connected to the 2nd, 3rd, 4th and so on, the second layer is connected to the 3rd, 4th, 5th and so on. This is done to enable maximum information flow between the layers of the network. To preserve the feed-forward nature, each layer obtains inputs from all the previous layers and passes on its own feature maps to all the layers which will come after it. Unlike Resnets it does not combine features through summation but combines the features by concatenating them. So the 'ith' layer has 'i' inputs and consists of feature maps of all its preceding convolutional blocks. Its own feature maps are passed on to all the next 'I-i' layers.

This introduces (I(I+1))/2 connections in the network, rather than just I' connections as in traditional deep learning architectures. It hence requires fewer parameters than traditional convolutional neural networks, as there is no need to learn unimportant feature maps.

DenseNet consists of two important blocks other than the basic convolutional and pooling layers. They are the Dense Blocks and the Transition layers.



DenseNet121 Model Building:

```
base_model = DenseNet121(include_top=False, weights="imagenet", input_shape=(224,224,3), pooling="avg")
```

model3 = Sequential()
model3.add(base_model)
model3.add(layers.Flatten())
model3.add(layers.Dense(2048 ,activation='relu'))
model3.add(BatchNormalization())

model3.add(Dropout(0.5))
model3.add(layers.Dense(2, activation ='sigmoid'))

model3.compile(optimizer='adam',loss='binary_crossentropy',metrics=['accuracy','Precision','Rec all'])

model3.summary()

Model: "sequential_1"

Layer (type)	Output Shape	Param #		
densenet121 (Function	onal) (None, 1024)	7037504		
flatten_2 (Flatten)	(None, 1024)	0		
dense_3 (Dense)	(None, 2048)	2099200		
batch_normalization (BatchN (None, 2048) 8192 ormalization)				
dropout_1 (Dropout)	(None, 2048)	0		
dense_4 (Dense)	(None, 2)	4098		

Total params: 9,148,994
Trainable params: 9,061,250
Non-trainable params: 87,744

history3 =


model3.fit(train_generator,steps_per_epoch=25,epochs=10,validation_data=val_generator)

```
Epoch 1/10
precision: 0.6990 - recall: 0.7025 - val_loss: 5.9835 - val_accuracy: 0.2290 - val_precision: 0.2290 -
val_recall: 0.2290
Epoch 2/10
25/25 [=============] - 222s 9s/step - loss: 0.7363 - accuracy: 0.7100 -
precision: 0.6897 - recall: 0.6725 - val_loss: 2.0226 - val_accuracy: 0.7690 - val_precision: 0.7353 -
val_recall: 0.0250
Epoch 3/10
precision: 0.6430 - recall: 0.6800 - val_loss: 26.1229 - val_accuracy: 0.2700 - val_precision: 0.2510 -
val_recall: 0.2440
Epoch 4/10
precision: 0.6574 - recall: 0.6525 - val_loss: 4.0967 - val_accuracy: 0.6550 - val_precision: 0.6574 -
val recall: 0.6620
Epoch 5/10
precision: 0.6863 - recall: 0.7000 - val_loss: 68.0669 - val_accuracy: 0.7580 - val_precision: 0.7608 -
val_recall: 0.7410
Epoch 6/10
precision: 0.7080 - recall: 0.6850 - val_loss: 9.5813 - val_accuracy: 0.4760 - val_precision: 0.4779 -
val recall: 0.4650
Epoch 7/10
precision: 0.7146 - recall: 0.7450 - val_loss: 1.0643 - val_accuracy: 0.5310 - val_precision: 0.5313 -
val_recall: 0.5350
Epoch 8/10
precision: 0.7410 - recall: 0.7225 - val_loss: 2.5277 - val_accuracy: 0.7720 - val_precision: 0.7718 -
val_recall: 0.7710
Epoch 9/10
precision: 0.7291 - recall: 0.7400 - val_loss: 0.7689 - val_accuracy: 0.7840 - val_precision: 0.7850 -
val_recall: 0.7850
Epoch 10/10
```

#Training Accuracy of the DenseNet121 Model score3 = getScore(history3)

The Average Training Accuracy is 0.7097499966621399

#Plots-accuracy and loss in train vs validation data plotAccAndLoss(history3)



#Test Data Metrics-Accuracy,Precision and Recall
result3 = model3.evaluate_generator(test_generator, 5000)
printTable([result3], headDataMetrics)

display table printTable(accData, headAccTable)

Model-1 VGG-16	0.718805	0.754	1
Model-2 ResNet50	0.6725	0.8024	1
Model-3 DenseNet121	0.70975	0.573	

Object Detection:

Object detection is a computer technology related to computer vision and image processing that deals with detecting instances of semantic objects of a certain class (such as humans, buildings, or cars) in digital images and videos. Well-researched domains of object detection include face detection and pedestrian detection. Object detection has applications in many areas of computer vision, including image retrieval and video surveillance.

For object detection, we have implemented a specific model whose steps are very clearly explained in the comments after each piece of code. First we load the needed packages and then import the dataset. Then, we split the data into train and validation samples. Next we define some functions and classes to be implemented in the model. This is tailored for the dataset on hand. Intersection Over Union metric and loss function are defined. The IOU score is the most important metric in object detection. It is the area of overlap divided by the area of the union. Generally, an IOU score of 0.5 and more is considered good for an object detection model, a threshold which is achieved by our model in the validation. After the model is built, we create a network and compile it with Adam optimizer, loss as iou_bce_loss as defined above and two metrics: accuracy and IOU score. Cosine function is used as the learning rate annealing function and 'generator' is used to create train and validation generators and later test generators. Finally, the model is fit with the parameters. Training vs validation metrics plots are visualized. At the end, the test data is used and the submission CSV file with predictions is given as an output. For further improvement, we can try different models like YOLO, Mask R-CNN, etc...

Object Detection Model Building:

The table contains [filename : pneumonia location] pairs per row.

A filename that has no pneumonia does not have a pneumonia location.

Likewise,a filename that contains multiple pneumonia contains multiple rows with the same filename but different pneumonia locations.

The table is loaded and transformed it into a dictionary.

The dictionary uses the filename as key and a list of pneumonia locations in that filename as value.

```
pneumonia_locations = {}
with open(os.path.join('stage_2_train_labels.csv'), mode='r') as infile:
  # open reader
  reader = csv.reader(infile)
  # skip header
  next(reader, None)
  # loop through rows
  for rows in reader:
    # retrieve information
    filename = rows[0]
    location = rows[1:5]
    pneumonia = rows[5]
    if pneumonia == '1':
      # convert string to float to int
      location = [int(float(i)) for i in location]
      # save pneumonia location in dictionary
      if filename in pneumonia_locations:
         pneumonia_locations[filename].append(location)
      else:
         pneumonia_locations[filename] = [location]
#load and shuffle filenames
folder = 'train_images'
filenames = os.listdir(folder)
random.shuffle(filenames)
#split into train and validation filenames
```

```
n_valid_samples = 2500
train_filenames = filenames[n_valid_samples:]
valid_filenames = filenames[:n_valid_samples]
print('train samples =', len(train_filenames))
print('validation samples =', len(valid_filenames))
n_train_samples = len(filenames) - n_valid_samples
train samples = 24218
validation samples = 2500
class generator(keras.utils.Sequence):
  def __init__(self, folder, filenames, pneumonia_locations=None, batch_size=32,
image_size=224, shuffle=True, augment=False, predict=False):
    self.folder = folder
    self.filenames = filenames
    self.pneumonia_locations = pneumonia_locations
    self.batch_size = batch_size
    self.image_size = image_size
    self.shuffle = shuffle
    self.augment = augment
    self.predict = predict
    self.on_epoch_end()
  def __load__(self, filename):
    # load dicom file as numpy array
    img = dicom.dcmread(os.path.join(self.folder, filename)).pixel_array
    # create empty mask
    msk = np.zeros(img.shape)
    # get filename without extension
    filename = filename.split('.')[0]
    # if image contains pneumonia
    if filename in pneumonia_locations:
```

```
# loop through pneumonia
    for location in pneumonia_locations[filename]:
       # add 1's at the location of the pneumonia
      x, y, w, h = location
      msk[y:y+h, x:x+w] = 1
  # if augment then horizontal flip half the time
  if self.augment and random.random() > 0.5:
    img = np.fliplr(img)
    msk = np.fliplr(msk)
  # resize both image and mask
  img = resize(img, (self.image_size, self.image_size), mode='reflect')
  msk = resize(msk, (self.image_size, self.image_size), mode='reflect') > 0.5
  # add trailing channel dimension
  img = np.expand_dims(img, -1)
  msk = np.expand_dims(msk, -1)
  return img, msk
def __loadpredict__(self, filename):
  # load dicom file as numpy array
  img = dicom.dcmread(os.path.join(self.folder, filename)).pixel_array
  # resize image
  img = resize(img, (self.image_size, self.image_size), mode='reflect')
  # add trailing channel dimension
  img = np.expand_dims(img, -1)
  return img
def __getitem__(self, index):
  # select batch
  filenames = self.filenames[index*self.batch_size:(index+1)*self.batch_size]
  # predict mode: return images and filenames
  if self.predict:
    # load files
```

```
imgs = [self.__loadpredict__(filename) for filename in filenames]
      # create numpy batch
      imgs = np.array(imgs)
      return imgs, filenames
    # train mode: return images and masks
    else:
      # load files
      items = [self.__load__(filename) for filename in filenames]
      # unzip images and masks
      imgs, msks = zip(*items)
      # create numpy batch
      imgs = np.array(imgs)
      msks = np.array(msks)
      return imgs, msks
  def on_epoch_end(self):
    if self.shuffle:
      random.shuffle(self.filenames)
  def __len__(self):
    if self.predict:
      # return everything
      return int(np.ceil(len(self.filenames) / self.batch_size))
    else:
      # return full batches only
      return int(len(self.filenames) / self.batch_size)
#Intersection Over Union or jaccard loss function
def iou_loss(y_true, y_pred):
  #print(y_true)
  y_true=tf.cast(y_true, tf.float32)
  y_pred=tf.cast(y_pred, tf.float32)
```

```
y_true = tf.reshape(y_true, [-1])
  y_pred = tf.reshape(y_pred, [-1])
  intersection = tf.reduce_sum(y_true * y_pred)
  score = (intersection + 1.) / (tf.reduce_sum(y_true) + tf.reduce_sum(y_pred) - intersection + 1.)
  return 1 - score
# combine bce loss and iou loss
def iou_bce_loss(y_true, y_pred):
  return 0.5 * keras.losses.binary_crossentropy(y_true, y_pred) + 0.5 * iou_loss(y_true, y_pred)
# mean iou as a metric
def mean_iou(y_true, y_pred):
  y_pred = tf.round(y_pred)
  intersect = tf.reduce_sum(y_true * y_pred, axis=[1, 2, 3])
  union = tf.reduce_sum(y_true, axis=[1, 2, 3]) + tf.reduce_sum(y_pred, axis=[1, 2, 3])
  smooth = tf.ones(tf.shape(intersect))
  return tf.reduce_mean((intersect + smooth) / (union - intersect + smooth))
def create_downsample(channels, inputs):
  x = keras.layers.BatchNormalization(momentum=0.9)(inputs)
  x = keras.layers.LeakyReLU(0)(x)
  x = keras.layers.Conv2D(channels, 1, padding='same', use_bias=False)(x)
  x = keras.layers.MaxPool2D(2)(x)
  return x
def create_resblock(channels, inputs):
  x = keras.layers.BatchNormalization(momentum=0.9)(inputs)
  x = keras.layers.LeakyReLU(0)(x)
  x = keras.layers.Conv2D(channels, 3, padding='same', use_bias=False)(x)
  x = keras.layers.BatchNormalization(momentum=0.9)(x)
  x = keras.layers.LeakyReLU(0)(x)
```

```
x = keras.layers.Conv2D(channels, 3, padding='same', use_bias=False)(x)
  return keras.layers.add([x, inputs])
def create_network(input_size, channels, n_blocks=2, depth=4):
  # input
  inputs = keras.Input(shape=(input_size, input_size, 1))
  x = keras.layers.Conv2D(channels, 3, padding='same', use_bias=False)(inputs)
  # residual blocks
  for d in range(depth):
    channels = channels * 2
    x = create_downsample(channels, x)
    for b in range(n_blocks):
      x = create_resblock(channels, x)
  # output
  x = keras.layers.BatchNormalization(momentum=0.9)(x)
  x = keras.layers.LeakyReLU(0)(x)
  x = keras.layers.Conv2D(128, 1, activation=None)(x)
  x = keras.layers.BatchNormalization(momentum=0.9)(x)
  x = keras.layers.LeakyReLU(0)(x)
  x = keras.layers.Conv2DTranspose(64, (8,8), (4,4), padding="same", activation=None)(x)
  x = keras.layers.BatchNormalization(momentum=0.9)(x)
  x = keras.layers.LeakyReLU(0)(x)
  x = keras.layers.Conv2D(1, 1, activation='sigmoid')(x)
  outputs = keras.layers.UpSampling2D(2**(depth-2))(x)
  model = keras.Model(inputs=inputs, outputs=outputs)
  return model
# create network and compiler
model = create_network(input_size=224, channels=32, n_blocks=2, depth=4)
model.compile(optimizer='adam', loss=iou_bce_loss, metrics=['accuracy', mean_iou])
```

```
# cosine learning rate annealing
def cosine_annealing(x):
 Ir = 0.001
 epochs = 7
 return lr*(np.cos(np.pi*x/epochs)+1.)/2
learning_rate = tf.keras.callbacks.LearningRateScheduler(cosine_annealing)
# create train and validation generators
folder = 'train_images'
train_gen = generator(folder, train_filenames, pneumonia_locations, batch_size=32,
image_size=224, shuffle=True, augment=False, predict=False)
valid_gen = generator(folder, valid_filenames, pneumonia_locations, batch_size=32,
image_size=224, shuffle=False, predict=False)
MULTI_PROCESSING = True
history = model.fit_generator(train_gen, validation_data=valid_gen, callbacks=[learning_rate],
use_multiprocessing=True, epochs=5,steps_per_epoch=100, shuffle=True, verbose=1)
Epoch 1/5
- mean iou: 0.6321 - val loss: 0.5386 - val accuracy: 0.9634 - val mean iou: 0.5955 - Ir: 0.0010
Epoch 2/5
100/100 [============================== ] - 5440s 54s/step - loss: 0.5207 - accuracy: 0.9543
- mean iou: 0.5617 - val loss: 0.5302 - val accuracy: 0.9254 - val mean iou: 0.5042 - Ir:
9.5048e-04
Epoch 3/5
- mean iou: 0.5810 - val loss: 0.4832 - val accuracy: 0.9635 - val mean iou: 0.6549 - Ir:
8.1174e-04
Epoch 4/5
- mean iou: 0.5971 - val loss: 0.4662 - val accuracy: 0.9637 - val mean iou: 0.6159 - Ir:
6.1126e-04
```

Epoch 5/5


```
100/100 [===============] - 5723s 57s/step - loss: 0.4636 - accuracy: 0.9621 - mean_iou: 0.6003 - val_loss: 0.4648 - val_accuracy: 0.9703 - val_mean_iou: 0.6755 - lr: 3.8874e-04
```

```
#Plots:Training vs validation-Loss,Accuracy and Intersection Over Union Score plt.subplot(131)
plt.plot(history.epoch, history.history["loss"], label="Train loss")
plt.plot(history.epoch, history.history["val_loss"], label="Valid loss")
plt.legend()

plt.subplot(132)
plt.plot(history.epoch, history.history["accuracy"], label="Train accuracy")
plt.plot(history.epoch, history.history["val_accuracy"], label="Valid accuracy")
plt.legend()

plt.subplot(133)
plt.plot(history.epoch, history.history["mean_iou"], label="Train iou")
```

plt.plot(history.epoch, history.history["mean_iou"], label="Train iou")
plt.plot(history.epoch, history.history["val_mean_iou"], label="Valid iou")
plt.legend()
plt.show()



load and shuffle filenames
folder = 'test_images'
test_filenames = os.listdir(folder)
print('test samples =', len(test_filenames))

```
# create test generator with predict flag set to True
test_gen = generator(folder, test_filenames, None, batch_size=32, image_size=224,
shuffle=False, predict=True)
# create submission dictionary
submission_dict = {}
# loop through testset
for imgs, filenames in tqdm(test_gen):
  # predict batch of images
  preds = model.predict(imgs)
  # loop through batch
  for pred, filename in zip(preds, filenames):
    # resize predicted mask
    pred = resize(pred, (1024, 1024), mode='reflect')
    # threshold predicted mask
    comp = pred[:, :, 0] > 0.5
    # apply connected components
    comp = measure.label(comp)
    # apply bounding boxes
    predictionString = "
    for region in measure.regionprops(comp):
      # retrieve x, y, height and width
      y, x, y2, x2 = region.bbox
      height = y2 - y
      width = x2 - x
      # proxy for confidence score
      conf = np.mean(pred[y:y+height, x:x+width])
      # add to predictionString
      predictionString += str(conf) + '' + str(x) + '' + str(y) + '' + str(width) + '' + str(height) + ''
    # add filename and predictionString to dictionary
    filename = filename.split('.')[0]
    submission_dict[filename] = predictionString
```

```
# stop if we've got them all
if len(submission_dict) >= len(test_filenames):
    break

print("Prediction complete")

# save dictionary as csv file-submission
sub = pd.DataFrame.from_dict(submission_dict,orient='index')
sub.index.names = ['patientId']
sub.columns = ['PredictionString']
sub.to_csv('submission.csv')

test samples = 3000
```

93/94 [27:11<00:17, 17.54s/it]

Prediction complete